

# AeroTrak® Cleanroom Condensation Particle Counter

**Model 9001** 



The AeroTrak® 9001 Cleanroom Condensation Particle Counter (CPC) is the only water-based, laminar flow, high flow rate instrument on the market—providing the confidence needed for monitoring critical environments.

Using patented technology, the AeroTrak® 9001 Cleanroom CPC provides 10 nm sensitivity at a 0.1 CFM (2.83 L/min) flow rate. With an ultra-low false count rate, the AeroTrak 9001 Cleanroom CPC is ideally suited for the low particle concentrations found in ISO Class 1 and Class 2 environments.

TSI is the world leader in condensation particle counter technology, supplying commercial instruments since 1979. TSI's world-class classifiers and CPCs are used by national and standards labs worldwide. This same technology is used in the AeroTrak 9001 Cleanroom CPC. Backed by a standard one-year warranty and TSI's long-standing reputation for high quality, there is no other CPC like it on the market today that provides the same confidence for monitoring critical environments.

Unlike other high flow rate water-based CPCs, the AeroTrak 9001 Cleanroom CPC does not use turbulent mixing, so it minimizes particle losses and increases counting accuracy. Unlike alcohol or other non-water based CPCs, the AeroTrak 9001 Cleanroom CPC does not introduce chemicals into clean environments.

The AeroTrak 9001 Cleanroom CPC has the same easy-to-use graphics interface as other AeroTrak Particle Counters. Multiple particle counting configurations can be conveniently stored and uploaded as needed.

The AeroTrak 9001 Cleanroom CPC can be used as a standalone particle counter or run continuously, integrated into a facility monitoring system like the TSI FMS Software. Optionally, the instrument can be interfaced to other systems using the Ethernet interface and Modbus communications.

#### **Features and Benefits**

- 10 nm sensitivity
- 0.1 CFM (2.83 l/min) sample flow
- Low false-count rate (< 0.07 count/ft³) provides accurate measurements at low particle concentrations
- Fast response
- Uses Ultra-Pure or DI water as working fluid, keeping alcohol or other chemicals out of the cleanroom
- Data Storage: 250 zones, 999 locations, and 10,000 samples
- USB, USB storage device, and Ethernet data transfer
- External reporting via TrakPro™ Lite Secure Software
- Industry's loudest audible alarm
- Easy, icon-driven touch screen operation
- Allows storage and recall of specific recipes and zone configurations
- Works with optional AeroTrak® 7950 High Pressure Diffuser (HPD) for measurement of high pressure gases (CDA, N<sub>2</sub>, AR)
- 316L stainless steel enclosure

## TSI—Providing the confidence needed for monitoring critical environments

- Ultra-low false count rate
- Trusted performance
  - Three detector design ensures reliability
- Comprehensive self-diagnostics



### **Specifications**

## AeroTrak® Condensation Particle Counter

		9001		
Particle Size Range	Min. Detectable Particle (D <sub>50</sub> ) - 10 nm		Unit ID	Configurable IP add
	Max. Detectable Particle - 1000 nm		Security	2-level password p
Particle Channel Sizes	Single channel			250 zones
Counting Efficiency	50% at 10 nm		999 Locations	
Concentration Limits	57,000,000 p/Cu. ft @ 10% coincidence		Data Storage	10,000 sample rec
Working Fluid	Ultra-Pure or De-Ionized Water		particle concentrationstrument status	
Light Source	Long life laser diode		Via USB storage de	
Zero Count	<0.07 count/Cu. ft		Data Transfer	Lite Secure Softwa cable, or optional T
Flow Rate	0.1 CFM (2.83 L/min)		Calibration Frequency	Recommended min
Sample Tube Extension	Up to 3 ft (1 m) for 10 nm particles	Dimension (H x W x Weight Power External Surface	Dimension (H x W x D)	22.4 in. x 9.3 in. x 8 (56.9 cm x 24.6 cm
Sampling Modes	Manual, automatic, beep; count or concentration		Woight	(without water bott 28.5 lbs (13.0 kg)
Sampling Time	1 second to 99 hours		· ·	, 0,
Sampling Frequency	1 to 9999 cycles or continuous			110 to 240 VAC ur
Vacuum Source	Internal pump		External Surface	316L Stainless ste
Communication Mode	Modbus TCP over Ethernet or USB		Operating Range	59° to 95° F (15° to 0 to 90% RH, non-
Audible Alarm	Built-in; > 85 dBA at 1 m (adjustable)		Storage Range	-40° to 158° F (-40
Alarm Output	Dry contacts, closed when alarm is engaged		Warranty	One year
Alarm Limits	Programmable		Included Accessories	Quick Start Guide, purge filter, 2 inlets
Display	VGA 5.7 in. (14.5 cm) touch screen display			OD straight tube co bracket for water b USB cable
Languages	English, German, French, Spanish, Japanese, Chinese (simplified), Italian			
Software	Compatible with TrakPro Lite Secure and FMS 5 Software		Optional Accessories	AeroTrak® 7950 Hi
				An assessment, increpresentative mod

01				
Unit ID	Configurable IP address			
Security	2-level password protection to lock out usage and configuration			
	250 zones			
Data Storage	999 Locations			
Data Storage	10,000 sample records including date, time, particle concentration, flow status and instrument status			
Data Transfer	Via USB storage device, connection to TrakPro Lite Secure Software through Ethernet of USB cable, or optional TSI FMS Software			
Calibration Frequency	Recommended minimum of once per year			
Dimension (H x W x D)	22.4 in. x 9.3 in. x 8.6 in. (56.9 cm x 24.6 cm x 21.8 cm) (without water bottle)			
Weight	28.5 lbs (13.0 kg)			
Power	110 to 240 VAC universal power supply			
External Surface	316L Stainless steel			
Operating Range	59° to 95° F (15° to 35° C) 0 to 90% RH, non-condensing			
Storage Range	-40° to 158° F (-40° to 70° C)			
Warranty	One year			
Included Accessories	Quick Start Guide, manual on a CD, power supply, purge filter, 2 inlets (barbed inlet and a 3/16 in. OD straight tube connector inlet), water bottle, bracket for water bottle, drain bottle and USB cable			
Optional Accessories	AeroTrak® 7950 High Pressure Diffuser			
ISO 14644-14:2016	An assessment, including an agreed upon representative mode of operation, according to ISO 14644-14 showed that the AeroTrak® 9001 Cleanroom CPC has suitability for use within a cleanroom of ISO Class 3 at 0.1µm			

Specifications are subject to change without notice.

AeroTrak, TSI and the TSI logo are registered trademarks of TSI Incorporated in the United States and may be protected under other country's trademark registrations.

To see our list of patents please visit: www.tsi.com/patents



**TSI Incorporated** - Visit our website **www.tsi.com** for more information.

Tel: +1 800 874 2811 Tel: +44 149 4 459200 Tel: +33 1 41 19 21 99 USA India UK China France Tel: +49 241 523030 Germany

Tel: +91 80 67877200 Tel: +86 10 8219 7688 **Singapore** Tel: +65 6595 6388

